



ARROWgard Blue® Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

REF **AK-22142-F**

14 Fr. | **2** Lumen | **15** cm catheter length | **.038** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 14 Fr. x 15 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .038" (0.97 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Pressure Transduction Probe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
- 1: Tissue Dilator: 16 Fr. (5.3 mm) x 11.4 cm
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®² One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 2: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Needle Holder
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Silk with Curved Needle

¹Licensed under US Patent No. 6,706,024.
²A registered trademark of CareFusion or one of its subsidiaries.



Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.


California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.

The ARROWgard Blue antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (10 Ga.)	1.40	250 400	70 100
Proximal (10 Ga.)	1.20	250 400	70 150

* Priming volumes are approximate and are done without accessories.
† Venous Pressure settings are approximate and were determined by varying the blood flow rates from 250 to 400 mL/min in an *in vitro* experiment using blood analog solution.

DuPont™ Tyvek®
made in USA

Manufacturer:
Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

IK-22142-100H, Rev. 01 (2020-12)



ARROWgard Blue® Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

REF **AK-22142-F**

14 Fr. | **2** Lumen | **15** cm catheter length | **.038** inch dia. spring-wire guide

