ARR()

Central Venous Catheterization Kit with Blue FlexTip[®] ARROWg⁺ard Blue[®] Catheter

16^{Ga.} 20^{cm}_{catheter} 0

1:

1:

1:

1:

1:

1:

1:

1:

1:

3:

1:

REF AK-24301

Contents:

CVC

- Indwelling Catheter: 16 Ga. x 20 cm Radiopaque Polyurethane with Blue 1: FlexTip[®], ARROWgard[®] Antimicrobial Surface Treatment¹, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft 1: Tip on One End - "J" Tip on Other) with Arrow Advancer
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer 1: Needle
- Injection Needle: 25 Ga. x 1" (2.54 cm) 1:
- 1: Local Anesthetic Components
- Pressure Transduction Probe 1:
- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson 1: Spring-Wire Introduction Syringe
- Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe 1:
- Syringe: 3 mL Luer-Lock 1:
- Cautior Consult

instructions for use

Do not reuse

STERILE EO Sterilized by ethylene oxide





Tissue Dilator

Scalpel: #11

5 mL 1% Lidocaine HCl Solution

SecondSite[™] Adjustable Hub: Fastener

Gauze Pad: 4" x 4" (10 cm x 10 cm)

SecondSite[™] Adjustable Hub: Catheter Clamp

Suture: 3-0 Braided Silk with Straight Needle

²A registered trademark of CareFusion or one of its subsidiaries.

Drape: 21" x 36" with 3" fenestration, with adhesive

with Hi-Lite Orange[™] Tint

SharpsAway® Disposal Cup

¹Licensed under US Patent No. 6,706,024.

Dust Cap: Non-Vented

Do not use if package is damaged

Not made with natural rubber latex

32 inch dia. spring-wire

3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution

Store between

20 - 25°C (68 - 77°F).

DuPont[™] Tyvek[®] | 🖉


Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Do not resterilize

All contents sterile unless otherwise indicated. Fluid path components are non-pyrogenic.

Contraindications: The ARROWg ard Blue® antimicrobial catheter technology is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



inch dia.

Manufacturer: Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



www.teleflex.com

IK-24301-1000 (3/17)



Central Venous Catheterization Kit with Blue FlexTip[®] ARROWg⁺ard Blue[®] Catheter



6^{Ga.} 20^{cm}_{catheter} .03 spring-wire

