



## Central Venous Catheterization Kit with Blue FlexTip® ARROWgard Blue® Catheter

REF AK-24306

**16** Ga. | **16** cm catheter length | **.025** inch dia. spring-wire guide

**Contents:**

- 1: Indwelling Catheter: 16 Ga. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment<sup>1</sup>, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Local Anesthetic Components
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator

- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 21" x 36" with 3" fenestration, with adhesive
- 1: Scalpel: #11
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>Licensed under US Patent No. 6,706,024.  
<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

DuPont™ Tyvek®  
HDPE  
packaging




Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

All contents sterile unless otherwise indicated. Fluid path components are non-pyrogenic.

**Contraindications: The ARROWgard Blue® antimicrobial catheter technology is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr.)
Distal (16 Ga.)	0.42	3550

\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

**Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

IK-24306-100K (3/17)



## Central Venous Catheterization Kit with Blue FlexTip® ARROWgard Blue® Catheter

REF AK-24306

**16** Ga. | **16** cm catheter length | **.025** inch dia. spring-wire guide