







Large-Bore Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip®, Arrowg+ard Blue® Catheter for **High Volume Infusions** 

cm catheter

## **Contents:**

- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One
- End "L" Tip on Other) with Arrow Advancer Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- Injection Needle: 25 Ga. x 1" (2.54 cm)
- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- Pressure Transduction Probe
- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow $^{\circ}$  Raulerson Spring-Wire Introduction Syringe Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- Syringe: 3 mL Luer-Lock
- Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep ®2 One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- SharpsAway® Disposal Cup
- Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- Gauze Pad: 4" x 4" (10 cm x 10 cm)
- Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>Licensed under US Patent No. 6,706,024. <sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

**⚠** California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfadrugs (refer to product instructions for references).

12 Fr. 40 mm	12 Ga.	7	
4.0 mm	12Ga.	J	

Lumen	(mL)	(mL/min)	
Distal (16 Ga.)	0.40	49	
Medial (12 Ga.)	1.30	196	
Proximal (12 Ga.)	1.20	221	
* Priming volumes are approximate.			

Priming Volume\*

Flow Rate<sup>†</sup>

imate. † Based on 100 cm head height. Catheter Cross Section





















Manufacturer: Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IK-25123-115A (2020-03)

**REF AK-25123-F** 







