



## Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® ARROWgard Blue® Catheter

REF AK-25703-A

**7** Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: 3 mL Luer-Lock Syringe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Tissue Dilator
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" with 4" fenestration
- 1: Local Anesthetic Components
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>Licensed under US Patent No. 6,706,024.  
<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

DuPont™ Tyvek® HDPE packaging


Caution	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Keep away from sunlight	Keep dry	Do not use if package is damaged	Not made with natural rubber latex	Store between 20 - 25°C (68 - 77°F).

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

All contents sterile unless otherwise indicated. Fluid path components are non-pyrogenic.

**Contraindications: The ARROWgard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr.)
Distal (16 Ga.)	0.40	2720
Medial (18 Ga.)	0.40	1260
Proximal (18 Ga.)	0.40	1390

\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

**Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**

www.teleflex.com

IK-25703-100Q (1/17)



## Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® ARROWgard Blue® Catheter

REF AK-25703-A

**7** Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide