



Quad-Lumen Central Venous Catheterization Kit with Blue FlexTip® ARROWg⁺ard Blue® Catheter

REF **AK-25854** **8.5** Fr. | **4** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- | | |
|--|--|
| <ul style="list-style-type: none"> 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard® Antimicrobial Surface Treatment, Extension Line Clamps 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™ 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle 1: Injection Needle: 25 Ga. x 1" (2.54 cm) 1: Syringe: 3 mL Luer-Lock 1: Pressure Transduction Probe 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe 1: Local Anesthetic Components | <ul style="list-style-type: none"> 1: Tissue Dilator 1: 5 mL 1% Lidocaine HCl Solution 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution with Hi-Lite Orange™ Tint 4: Dust Cap: Non-Vented 1: SecondSite™ Adjustable Hub: Fastener 1: SecondSite™ Adjustable Hub: Catheter Clamp 1: SharpsAway® Disposal Cup 1: Drape: 24" x 36" with 4" fenestration 1: Scalpel: #11 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) 1: Suture: 3-0 Braided Silk with Straight Needle |
|--|--|

¹A registered trademark of CareFusion or one of its subsidiaries.

DuPont™ Tyvek®
packaging

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.


All contents sterile unless otherwise indicated.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Fluid path components are non-pyrogenic.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter technology is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Caution	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Keep away from sunlight	Keep dry	Do not use if package is damaged	Not made with natural rubber latex	Store between 20 - 25°C (68 - 77°F).



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr.)
Distal (16 Ga.)	0.43	2700
Medial 1 (14 Ga.)	0.54	5100
Medial 2 (18 Ga.)	0.36	1400
Proximal (18 Ga.)	0.41	1600

* Priming volumes are approximate and are done without accessories.
† Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IK-25854-100Q (2/17)



Quad-Lumen Central Venous Catheterization Kit with Blue FlexTip® ARROWg⁺ard Blue® Catheter

REF **AK-25854** **8.5** Fr. | **4** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide