

REF AK-29803-CDC

ARROW



Percutaneous Sheath Introducer Kit with Arrowgard Blue® Sheath and Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

8.5 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Arrowgard® Antimicrobial Surface Treatment<sup>1</sup>, Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: Dust Cap: Non-Vented
- 1: SharpsAway II™ Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Drape: 41" x 55" with 4" fenestration

- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Suture: 3-0 Braided Silk with Straight Needle
- 1: Tubing: 20" (50 cm) Extension

<sup>1</sup>Licensed under US Patent Nos. 6,706,024.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.


⚠ **California Prop. 65**

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.

**Contraindications:** The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

**Warning:** The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.



Not made with natural rubber latex

<sup>≠</sup>Precaution: Contains Phthalate: DEHP

Store between 20-25°C (68-77°F).

Manufacturer: Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA

Teleflex® [www.teleflex.com](http://www.teleflex.com)

IK-29803-105A (12/19)

REF AK-29803-CDC

ARROW



Percutaneous Sheath Introducer Kit with Arrowgard Blue® Sheath and Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

8.5 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide