



ARROWg⁺ard Blue PLUS[®] Quad-Lumen CVC

REF **AK-42854-CDC**

8.5 Fr. | **4** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip[®], ARROWg⁺ard Blue Plus[®] Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Components
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock

- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep^{®2} One-Step Solution with Hi-Lite Orange[™] Tint
- 4: Dust Cap: Non-Vented
- 1: SecondSite[™] Adjustable Hub: Fastener
- 1: SecondSite[™] Adjustable Hub: Catheter Clamp
- 1: SharpsAway II[™] Locking Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Safety Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.
²A registered trademark of CareFusion or one of its subsidiaries.

DuPont™ Tyvek™
packaging

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.


California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: The ARROWg⁺ard Blue PLUS[®] antimicrobial catheter technology is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

| | | | | | | | | | |
|---------|------------------------------|--------------|--------------------|------------------------------|-------------------------|----------|----------------------------------|------------------------------------|--------------------------------------|
| | | | | | | | | | |
| Caution | Consult instructions for use | Do not reuse | Do not resterilize | Sterilized by ethylene oxide | Keep away from sunlight | Keep dry | Do not use if package is damaged | Not made with natural rubber latex | Store between 20 - 25°C (68 - 77°F). |



| Lumen | Priming Volume* (mL) | Gravity Flow Rate [†] (mL/hr.) |
|-------------------|----------------------|---|
| Distal (16 Ga.) | 0.39 | 3200 |
| Medial 1 (14 Ga.) | 0.47 | 5400 |
| Medial 2 (18 Ga.) | 0.36 | 1600 |
| Proximal (18 Ga.) | 0.38 | 1800 |

* Priming volumes are approximate and are done without accessories.
[†] Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex[®]
www.teleflex.com

IK-42854-101F (5/17)



ARROWg⁺ard Blue PLUS[®] Quad-Lumen CVC



REF **AK-42854-CDC**

8.5 Fr. | **4** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide