



ARROWgard Blue Plus[®] Multi-Lumen CVC Kit

REF AK-45703-A

7 Fr. | 3 Lumen | 20 cm catheter length | .032 inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip[®], ARROWgard Blue PLUS[®] Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer[™]
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Local Anesthetic Components
- 1: Tissue Dilator
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep^{®2} One-Step Solution with Hi-Lite Orange[™] Tint
- 3: Dust Cap: Non-Vented
- 1: SecondSite[™] Adjustable Hub: Fastener
- 1: SecondSite[™] Adjustable Hub: Catheter Clamp
- 1: SharpsAway[®] Disposal Cup
- 1: Drape: 24" x 36" with 4" fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2"
- 5: Gauze Pad: 4" x 4"
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.
²A registered trademark of CareFusion or one of its subsidiaries.

DuPont™ Tyvek[®] packaging

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Fluid path components are non-pyrogenic.

Contraindications: The ARROWgard Blue PLUS[®] antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, chlorhexidine acetate, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry




Do not use if package is damaged



Not made with natural rubber latex



Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F)



Lumen	Priming Volume*(mL)	Flow Rate† (mL/hr.)
Distal (16 Ga.)	0.44	3100
Medial (18 Ga.)	0.39	1500
Proximal (18 Ga.)	0.39	1600

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex[®]
 www.teleflex.com

IK-45703-101H (2/17)



ARROWgard Blue Plus[®] Multi-Lumen CVC Kit

REF AK-45703-A

7 Fr. | 3 Lumen | 20 cm catheter length | .032 inch dia. spring-wire guide