



ARROWgard Blue Plus® Multi-Lumen CVC Extreme Kit

REF **AK-45703-AX**

7 Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard Blue PLUS® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide™² 25 Ga. x 5/8" (1.60 cm)
- 1: 3 mL Luer-Lock Syringe
- 1: Injection Needle: SafetyGlide™² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Local Anesthetic Components
- 1: Tissue Dilator

- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®³ One-Step Solution with Hi-Lite Orange™ Tint
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Safety Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: STATLOCK®⁴ CV Plus Catheter Stabilization Device
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

²A trademark of Becton, Dickinson and Company.

³A registered trademark of CareFusion or one of its subsidiaries.

⁴A registered trademark of C. R. Bard, Inc.

DuPont™ Tyvek®
packaging

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Fluid path components are non-pyrogenic.

Contraindications: The ARROWgard Blue PLUS® antimicrobial catheter technology is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight




Keep dry




Do not use if package is damaged



Not made with natural rubber latex



Store between 20 - 25°C (68 - 77°F).



Catheter Cross Section

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr.)
Distal (16 Ga.)	0.44	3100
Medial (18 Ga.)	0.39	1500
Proximal (18 Ga.)	0.39	1600

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IK-45703-109L (1/17)



ARROWgard Blue Plus® Multi-Lumen CVC Extreme Kit

REF **AK-45703-AX**

7 Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide