







## Arrowg+ard Blue PLUS® Multi-Lumen CVC





20 cm catheter length

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard Blue PLUS® Antimicrobial Surface Treatment¹, Extension Line Clamps

  1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End "J" Tip on
- Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle 1: Injection Needle: 25 Ga. x 1" (2.54 cm) 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

- Pressure Transduction Probe
  Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction
- Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- Syringe: 3 mL Luer-Lock
  Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 5 mL 1% Lidocaine HCl Solution and Alcohol Prep 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)

- (non-sterile solution)
  3: Dust Cap: Non-Vented
  1: SecondSite™ Adjustable Hub: Fastener
  1: SecondSite™ Adjustable Hub: Catheter Clamp
  1: SharpsAway II™ Locking Disposal Cup
  1: SharpsAway® Disposal Cup
  1: Drape: 41" x 55" with 4" fenestration
  1: Safety Scalpel: #11

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

📤 California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.


Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue PLUS antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)

  1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>Licensed under US Patent Nos. 6.706.024, 6.872,195 and 7.329.412. <sup>2</sup>A registered trademark of CareFusion or one of its subsid



Lumen	(mL)	(mL/hr)
Distal (16 Ga.)	0.44	3100
Medial (18 Ga.)	0.39	1500
Proximal (18 Ga.)	0.39	1600

Priming volumes are approximate and are done without accessories. Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Priming Gravity
Volume\* Flow Rate†





















Manufacturer Arrow International, Inc Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IK-45703-116H (12/19)

REF AK-45703-ACDC







