

REF ML-00701

Central Venous Catheter

with Blue FlexTip®

for use only with Arrow MAC™

Two-Lumen Central Venous Access Device

ARROW®

Catheter Diameter: 7 Fr., 2.4 mm

Priming Volume: 0.61 mL

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

Flow Rate: 9865 mL/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.


Contents:

1: Indwelling Catheter: 7 Fr.
Radiopaque Polyurethane with Blue FlexTip®,
Catheter Contamination Guard,
Extension Line Clamp, Injection Site Cap

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not re-sterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

 Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Athlone, Ireland



- FR** Cathéter veineux central avec cathéter Blue FlexTip à utiliser exclusivement avec le dispositif d'accès veineux central à deux lumières Arrow MAC
- DE** Zentraler Venenkatheter mit Blue FlexTip nur zum Gebrauch mit dem Einführbesteck für den zweilumigen, zentralen Venenkatheter von Arrow MAC
- IT** Catetere per vena centrale con Blue FlexTip da usare esclusivamente con il dispositivo di accesso venoso centrale a due lumi Arrow MAC
- PL** Cewnik do cewnikowania żył centralnych z końcówką Blue FlexTip, do użytku tylko z dwukanalowym urządzeniem Arrow MAC do uzyskiwania dostępu do żył centralnych
- PT** Cateter Venoso Central com Blue FlexTip para utilizar apenas com o Dispositivo para Acesso Venoso Central com Lúmen Duplo Arrow MAC
- RU** Катетер для главных вен с Blue FlexTip, для применения исключительно с двухсветовым устройством для доступа в главные вены Arrow MAC
- SL** Centralni veniski kateter z Blue FlexTip je samo za uporabo s pripomočkom za veniski dostop Arrow MAC z dvema svetlinama
- ES** Catéter para vena central con Blue FlexTip para uso exclusivo con el dispositivo de acceso venoso central de dos luces de Arrow MAC
- SV** Centralvenös kateter med Blue FlexTip, endast för användning med Arrow MAC tvåkanalig produkt för centralt ventillträde
- TR** Santral Venöz Kateter, Blue FlexTip ile, sadece Arrow MAC İki Lümenli Santral Venöz Erişim Cihazı ile kullanılmak içindir

IL-00701-1308 (2/12)

ARROW
INTERNATIONAL

REF
ML-00701

Central Venous Catheter with Blue FlexTip®

7Fr.

14 cm indwelling
catheter length
with **MAC**™

Peel to
Open