


Lumen	Priming Vol.* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.40	4824
Proximal (16 Ga.)	0.46	4494

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Contents:

1: Two-Lumen Indwelling Catheter:
7 Fr. Radiopaque Polyurethane
with Blue FlexTip®, Catheter Contamination Guard,
Extension Line Clamps, Injection Site Caps

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

 Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Athlone, Ireland



REF ML-00702

Two-Lumen Central Venous Catheter

with Blue FlexTip®
for use only with Arrow MAC™
Two-Lumen Central Venous Access Device

ARROW®

- FR** Cathéter veineux central à deux lumières avec Blue FlexTip à utiliser exclusivement avec le dispositif d'accès veineux central à deux lumières Arrow MAC
- DE** Zweilumiger, zentraler Venenkatheter mit Blue FlexTip nur zum Gebrauch mit dem Einführbesteck für den zweilumigen, zentralen Venenkatheter von Arrow MAC
- IT** Catetere per vena a due lumi con Blue FlexTip da usare esclusivamente con il dispositivo di accesso venoso centrale a due lumi Arrow MAC
- PL** Dwukanalowy cewnik do cewnikowania żył centralnych z końcówką Blue FlexTip, do użytku tylko z dwukanalowym urządzeniem Arrow MAC do uzyskiwania dostępu do żył centralnych
- PT** Cateter Venoso Central com Lúmen Duplo e Blue FlexTip para utilizar apenas com o Dispositivo para Acesso Venoso Central com Lúmen Duplo Arrow MAC
- RU** Двухпросветный катетер для главных вен с Blue FlexTip, для применения исключительно с двухпросветным устройством для доступа главных вен Arrow MAC
- SL** Centralni venski kateter z dvema svetlinama z Blue FlexTip je samo za uporabo s pripomočkom za venski dostop Arrow MAC z dvema svetlinama
- ES** Catéter venoso central de dos luces con punta Blue FlexTip para uso exclusivo con el dispositivo de acceso venoso central de dos luces de Arrow MAC
- SV** Tvåkanalig central venkateter med Blue FlexTip, endast för användning med Arrow MAC tvåkanalig produkt för centralt ventillträde
- TR** İki Lümenli Santral Venöz Kateter, Blue FlexTip ile, sadece Arrow MAC İki Lümenli Santral Venöz Erişim Cihazı ile kullanılmak içindir

IL-00702-1308 (10/11)

ARROW
INTERNATIONAL

REF
ML-00702

7 Fr. 2 Lumen

Two-Lumen Central Venous Catheter

16 cm indwelling
catheter length
with **MAC™**

Peel to
Open