

## Arterial Arterial Catheterization Kit

## REF ASK-01618-CHUM

## Contents:

- 1: Indwelling Catheter: 18 Ga. x 6-1/2" (16 cm) Radiopaque Polyurethane
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One
- End "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW
- 1: Injection Needle: SafetyGlide<sup>™1</sup> 18 Ga. x 1-1/2" (3.81 cm) RW
- 1: Injection Needle: SafetyGlide<sup>™1</sup> 25 Ga x 1" (2.54 cm)
- 1: Syringe: 10 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: Tubing: 24" Pressure Monitoring<sup>±</sup>
- 1: SharpsAway II<sup>™</sup> Locking Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 4: Gauze Pad: 4" x 4"
- 1: Suture: 3-0 Silk with Curved Needle

<sup>1</sup> A trademark of Becton, Dickinson and Company



Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. Sterile. Fluid path components are non-pyrogenic.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Contraindications:

Inadequate collateral arterial circulation.

## **Neleflex**®

Arrow International, Inc. Subsidiary of Teleflex, Incorporated 2400 Bernville Road Reading, PA 19605 USA

www.teleflex.com

IM-01618-100A (11/15)



Arterial Catheterization Kit

REF ASK-01618-CHUM 18<sup>Ga.</sup> 16<sup>Ga.</sup> 1025



