

REF ASK-04018-UPM

ARROW®



Arterial Line Catheterization Kit with Sharps Safety Features

18^{Ga.} | 16 cm catheter length | .025 inch dia. spring-wire guide

Contents:

- 1: Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 18 Ga. x 3" (7.62 cm) RW
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 2: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®² One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: SharpsAway II™ Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Needle Holder
- 1: Scissors
- 1: Forceps
- 1: Safety Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 2: Suture: 2-0 Silk with Curved Needle
- 1: Tubing±: 20" (50 cm) Extension

¹A trademark of Becton, Dickinson and Company.

²A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

⚠ **California Prop. 65**

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.



STERILE EO



Not made with natural rubber latex



±Precaution: Contains Phthalate: DEHP



Store between 20 - 25°C (68 - 77°F).

REF ASK-04018-UPM

ARROW®



Arterial Line Catheterization Kit with Sharps Safety Features

18^{Ga.} | 16 cm catheter length | .025 inch dia. spring-wire guide

Manufacturer:
Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

IM-04018-105A, Rev. 01 (2020-12)