



## **Femoral Arterial Line Catheterization Kit**



Contents

- 1:
- Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other) Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle 1:
- 1:
- 1: 1: 1:
- Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW Introducer Needle: 18 Ga. x 3" (7.62 cm) RW
- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Syringe: 5 mL Lue-Slip 5 mL 1% Lidocaine HCl Solution and Alcohol Prep 1:
- 2: 1:
- Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration Scalpel: #11 1:
- 1: 1: Prep Razor
- 2.
- 4:
- 8:
- Prep Razor Prep Sponge Swab Gauze Pad: 2" x 2" (5 cm x 5 cm) Gauze Pad: 4" x 4" (10 cm x 10 cm) Suture: 3-0 Braided Silk with Straight Needle 1:
- Tray: Prep 1:

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. △ California Prop. 65 Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Manufacturer: Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA

ARROW









18<sup>Ga.</sup> 16<sup>cm</sup><sub>catheter</sub> 025<sup>inch dia.</sup> guide