



# Arterial Catheterization Kit



REF ASK-04510-UPM

20<sup>Ga.</sup> | 12<sup>cm</sup> catheter length | .025<sup>inch dia.</sup> spring-wire guide

## Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Components
- 2: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: SharpsAway II™ Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Needle Holder
- 1: Scissors

- 1: Forceps
- 1: Safety Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 2: Suture: 2-0 Silk with Curved Needle
- 1: Tubing±: 20" (50 cm) Extension

<sup>1</sup> A trademark of Becton, Dickinson and Company.  
<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.



Rx only. **Warning:** Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65  
Warning: Cancer and Reproductive Harm [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

**Warning:** The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
Fluid path components are non-pyrogenic.

**Manufacturer:**  
 Arrow International LLC  
 Subsidiary of Teleflex, Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IM-04510-101G, Rev. 01 (2020-12)



# Arterial Catheterization Kit

REF ASK-04510-UPM

20<sup>Ga.</sup> | 12<sup>cm</sup> catheter length | .025<sup>inch dia.</sup> spring-wire guide

