



# Percutaneous Sheath Introducer Kit

for use with 7.5 - 8 Fr. Catheters



REF **ASK-09903-SM**

**9** Fr. | **10** cm sheath length | **.035** inch dia. spring-wire guide

**Contents:**

- |   |   |  |
|---|---|--|
| <ul style="list-style-type: none"> <li>1: Sheath: 9 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator</li> <li>1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer</li> <li>1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle</li> <li>1: Pressure Transduction Probe</li> <li>1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe</li> <li>1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe</li> <li>1: Injection Needle: SafetyGlide™ 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe</li> <li>1: Syringe: 10 mL Luer-Lock</li> <li>1: Obturator: 8 Fr.</li> </ul> | <ul style="list-style-type: none"> <li>1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter</li> <li>1: Dust Cap: Non-Vented</li> <li>1: SharpsAway® Disposal Cup</li> <li>1: SharpsAway II™ Locking Disposal Cup</li> <li>1: Maximal Barrier Drape™ with 4" fenestration</li> <li>1: Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive</li> <li>1: Towel</li> <li>1: Needle Holder</li> <li>1: Forceps</li> <li>1: Safety Scalpel: #11</li> <li>1: Stopcock: 4-way High-Flow</li> <li>1: Checklist/CLIP Sheet</li> <li>1: Sterile Procedure Sign</li> <li>2: Gauze Pad: 2" x 2" (5 cm x 5 cm)</li> </ul> | <ul style="list-style-type: none"> <li>5: Gauze Pad: 4" x 4" (10 cm x 10 cm)</li> <li>1: Dressing: Tegaderm®<sup>2</sup> 10 cm x 12 cm</li> <li>1: Suture: 3-0 Braided Silk with Straight Needle</li> <li>1: Suture: 3-0 Silk with Curved Needle</li> <li>1: Tubing<sup>±</sup>: 20" (50 cm) Extension</li> <li>1: HemoHopper® Fluid Receptacle</li> </ul> |
|---|---|--|

<sup>1</sup> A trademark of Becton, Dickinson and Company.  
<sup>2</sup> A registered trademark of 3M Company.

**Warning:** The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.



Not made with natural rubber latex



±Precaution: Contains Phthalate: DEHP