



ARROW⁺gard Blue[®] PSI Kit for use with 7.5 - 8 Fr. Catheters



REF **ASK-29804-UHN2**

Contents:

- | | | |
|---|---|---|
| 1: Sheath: 8.5 Fr. x 4-1/8" (10 cm) Arrow-Flex [®] Radiopaque Polyurethane with ARROW ⁺ gard [®] Antimicrobial Surface Treatment ¹ , Integral Hemostasis Valve/Side Port and Tissue Dilator | 1: Injection Needle: SafetyGlide ^{™2} 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 1: Stopcock: 4-way High-Flow |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer [™] | 1: Syringe: 10 mL Luer-Lock | 1: Checklist/CLIP Sheet |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Obturator: 8 Fr. | 1: Sterile Procedure Sign |
| 1: Pressure Transduction Probe | 1: Cath-Gard [®] Catheter Contamination Shield: 80 cm with TwistLock [™] Adapter | 2: Gauze Pad: 2" x 2" |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow [®] Raulerson Spring-Wire Introduction Syringe | 1: Dust Cap: Non-Vented | 5: Gauze Pad: 4" x 4" |
| 1: Injection Needle: SafetyGlide ^{™2} 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe | 1: SharpsAway [®] Disposal Cup | 1: Dressing: Tegaderm ^{®3} 10 cm x 12 cm |
| | 1: SharpsAway II [™] Locking Disposal Cup | 1: Suture: 3-0 Silk with Curved Needle |
| | 1: Maximal Barrier Drape [™] with 4" fenestration | 1: Tubing: 20" Extension \pm |
| | 1: Drape: 21" x 36" with 3" fenestration, with adhesive | 1: HemoHopper [®] Fluid Receptacle |
| | 1: Towel | |
| | 1: Needle Holder | |
| | 1: Safety Scalpel: #11 | |

¹ Licensed under US Patent No. 6,706,024.

² A trademark of Becton, Dickinson and Company.

³ A registered trademark of 3M Company.

Contraindications: The ARROW⁺gard Blue[®] antimicrobial sheath is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

~~LATEX~~

Not made with
natural rubber latex

PHT DEHP

\pm Precaution:
Contains Phthalate: DEHP



Store below 25°C (77°F)
Avoid excessive heat above 40°C (104°F)