

Arterial Catheterization Kit



ARROW
INTERNATIONAL

Arrow Select Kits

ASK-00520-WHO

> **Sharps Safety Features**

Contents:

- 1: Indwelling Catheter: 20 Ga. x 2" (5 cm) Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
- 1: **Injection Needle: SafetyGlide^{®1} 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe**
- 1: Syringe: 5 mL Luer-Slip
- 1: **SharpsAway II[™] Locking Disposal Cup**
- 1: CSR Wrap
- 1: Drape: 30" x 36" with 3" fenestration, with adhesive
- 1: Dressing: Tegaderm^{®2} 10 cm x 12 cm
- 1: **Dressing: STATLOCK^{®3} PICC Plus Catheter Stabilization Device (prep pad solution non-sterile)**
- 3: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 1: Swabstick Foil Skin Protectant (non-sterile solution)

¹ A registered trademark of Becton, Dickinson and Company.

² A registered trademark of 3M Company.

³ A registered trademark of C. R. Bard, Inc.

Safety Features:

- > SafetyGlide[®] Protected Needle
- > SharpsAway II[™] Locking Disposal Cup
- > STATLOCK[®] Catheter Stabilization Device

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Contraindications: Inadequate collateral arterial circulation.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road, Reading, PA 19605 USA
www.teleflex.com

10-00520-100B (10/17)

Arterial Catheterization Kit

ASK-00520-WHO

20 Ga.

5 cm catheter length

.021 inch diameter spring-wire guide

> **Sharps Safety Features**

Arrow Select Kits



ARROW
INTERNATIONAL