

ARROWgtard Blue PLUS® Multi-Lumen CVC

REF MTO-42703-AM

Antimicrobial surfaces utilizing chlorhexidine, chlorhexidine acetate and silver sulfadiazine

Special Order Product

Multi-Lumen Indwelling Catheter: 7 Fr. x 6" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWgtard Blue PLUS® Antimicrobial Surface Treatment¹

Extension Line Clamps, Injection Site Caps

One: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow AdvancerTM

Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle

Pressure Transduction Probe One: Syringe: 5 mL Luer-Slip One:

Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe One:

Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe

One: Tissue Dilator

Three: UserGard® Needle-Free Injection Hubs

Clamp: Catheter Fastener: Catheter Clamp One:

One: SharpsAway IITM Locking Disposal Cup

Drape: 41" x 50" with 4" fenestration, with adhesive One:

One: Safety Scalpel: #11 Four: Staple Anchoring Devices^{2,3} Gauze Pads: 2" x 2' Two: Five: Gauze Pads: 4" x 4"


Suture: 3-0 Braided Silk with Straight Needle

INFECTION PROTECTION FEATURES:

- AGB^{+®} catheters now include internal lumen protection & enhanced external surface activity
- Extra-Large Drape size

SAFETY FEATURES:

- SharpsAway IITM
 Locking Disposal Cup
- Staple Anchoring Device Safety Scalpel



| Catheter Cross | Section |
|----------------|---------|

| Lumen | Priming Volume* (mL) | Flow Rate [†] (mL/hr) |
|-------------------|----------------------------|--------------------------------------|
| Distal (16 Ga.) | 0.38 | 3000 |
| Medial (18 Ga.) | 0.37 | 1400 |
| Proximal (18 Ga.) | 0.39 | 1600 |

^{*}Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL. flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Contraindications: The ARROWg[†]ard Blue PLUS® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, chlorhexidine acetate, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

³MR Conditional.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°













Arrow International, Inc. 2400 Bernville Road Reading, PA 19605 USA www.arrowintl.com







ARROWgtard Blue PLUS® Multi-Lumen CVC

spring-wire

¹Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412

²Licensed under Us Patent No. RE37,814.