



## PICC Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter Seldinger Access Set with 45 cm Wire

REF **PR-05041**

**4** Fr. | **50** cm catheter length | **.018** inch dia. spring-wire guide

**Contents:**

- 1: Indwelling Catheter: 4 Fr. x 50 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: GlideThru™ Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 1: Skin Protectant Prep Pad (non-sterile solution)
- 1: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure
- 1: Dressing: STATLOCK® Catheter Stabilization Device

\*A registered trademark of C. R. Bard, Inc.

All components are CE 0086 unless otherwise noted.



Not made with natural rubber latex

- fr** Dispositif d'accès Seldinger pour cathéter central inséré périphériquement (CCIP) avec cathéter Blue FlexTip avec fil de 45 cm
- de** PEZK-Besteck für peripher eingeführten, zentralen Katheter mit Blue FlexTip-Katheter für Seldinger-Zugang mit Draht (45 cm)
- it** Set di accesso Seldinger per catetere centrale a inserimento periferico, con catetere Blue FlexTip e filo da 45 cm
- pt** Conjunto de Cateter Central de Introdução Periférica PICC com Cateter Blue FlexTip por Técnica de Seldinger com Fio de 45 cm
- es** Conjunto de acceso Seldinger para catéter de vena central de inserción periférica (PICC) con catéter Blue FlexTip y alambre de 45 cm
- sv** Seldingerinföringsset för PICC perifert införd centralkateter med Blue FlexTip-kateter och 45 cm lång tråd


**Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.50	850	3700

\* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates determined at maximum pump pressure of 10 psig.

**Manufacturer:**

Arrow International, Inc.

Subsidiary of Teleflex Incorporated

2400 Bernville Road

Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

IR-05041-135A (04/19)

**CE** REP 0086 Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland



PICC Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter Seldinger Access Set with 45 cm Wire

REF **PR-05041**

**4** Fr. | **50** cm catheter length | **.018** inch dia. spring-wire guide