

# ARROW®



# ARROW® PICC Set

## REF PR-35052-HPSL

### Contents:

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>1: Two-Lumen Indwelling TaperFree™ Catheter: 5 Fr. x 50 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip® and Extension Line Clamps</li> <li>1: GlideThru™ Peel-Away Sheath: 5 Fr. x 4" (10 cm) Radiopaque over 5 Fr. Dilator</li> <li>1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)</li> <li>1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 51-3/16" (130 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)</li> <li>1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW</li> <li>1: Syringe: 5 mL Luer-Lock</li> <li>1: Skin Protectant Prep Pad</li> <li>2: Dust Caps: Non-Vented</li> </ul> | <ul style="list-style-type: none"> <li>1: SecondSite™ Adjustable Hub: Fastener</li> <li>1: SecondSite™ Adjustable Hub: Catheter Clamp</li> <li>1: Protection Sleeve</li> <li>1: Flow Rate Information Card with Injection Log</li> <li>1: Patient ID Card</li> <li>1: Chart Sticker</li> <li>1: Patient Information Booklet</li> <li>1: Checklist/CLIP Sheet</li> <li>1: Paper Tape Measure</li> <li>1: Dressing: STATLOCK®<sup>1</sup> Catheter Stabilization Device</li> </ul> |
|---|--|

<sup>1</sup>A registered trademark of C. R. Bard, Inc.  
All components are CE 0086 unless otherwise noted.



Caution



Consult  
instructions for use



Do not reuse



Do not resterilize



Sterilized  
by ethylene oxide



Keep away  
from sunlight



Keep dry



Do not use  
if package is damaged




Not made with  
natural rubber latex

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

**Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of patient must be completed to ensure no contraindications exist. See additional labeling for product specific contraindications.**



Manufacturer:  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com



IR-35052-104A, Rev. 2 (10/16)



# ARROW® PICC Set

## REF PR-35052-HPSL

5 Fr.

2 Lumen

50 cm  
catheter  
length

.018 inch dia.  
spring-wire  
guide

130 cm Stainless  
Steel Wire

# ARROW®