



## ARROW® JACC with Chloragard® Technology

### REF JR-42052-HPHNM

**Contents:**

- |   |   |
|---|---|
| <ul style="list-style-type: none"> <li>1: Two-Lumen Indwelling TaperFree™ CG+™ Antimicrobial/Antithrombogenic Catheter<sup>1</sup>: 5.5 Fr. x 20 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps and Contamination Guard</li> <li>1: GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator</li> <li>1: Guidewire, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)</li> <li>1: Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil (CE 0120)</li> <li>1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW</li> <li>1: Syringe: 10 mL Luer-Lock</li> <li>1: Skin Protectant Prep Pad</li> </ul> | <ul style="list-style-type: none"> <li>2: Dust Caps: Non-Vented</li> <li>1: SecondSite™ Adjustable Hub: Fastener</li> <li>1: SecondSite™ Adjustable Hub: Catheter Clamp</li> <li>1: Safety Scalpel: #11</li> <li>1: Patient ID Card</li> <li>1: Chart Sticker</li> <li>1: Patient Information Booklet</li> <li>1: Flow Rate Information Card with Injection Log</li> <li>1: Paper Tape Measure</li> <li>1: Dressing: STATLOCK® Catheter Stabilization Device</li> </ul> |
|---|---|

<sup>1</sup> Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

<sup>2</sup> A registered trademark of C. R. Bard, Inc.

**All components are CE 0086 unless otherwise noted.**



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex



Store below 25°C (77°F). Avoid excessive heat above 30°C (86°F).

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

**Contraindications: The Arrow® JACC with Chloragard® Technology is contraindicated for patients with known hypersensitivity to chlorhexidine (refer to product instructions for references).**



Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Manufacturer:  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com