







ARROW PICC with Chloragtard Technology

REF PR-45563-HPHNM

Contents:

- 1: Multi-Lumen Indwelling TaperFree™ CG+™
 Antimicrobial/Antithrombogenic Catheter¹: 6 Fr. x 55 cm
 Pressure Injectable Radiopaque Polyurethane with
 Blue FlexTip®, Extension Line Clamps
 and Contamination Guard
- 1: GlideThru™ Peel-Away Sheath: 6 Fr. x 4" (10 cm) Radiopague over 6 Fr. Dilator
- 1: Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil
- 1: Guidewire, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) Nitinol (Straight Soft Tip on One End – Straight Stiff Tip on Other)
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 10 mL Luer-Lock
- 1: Skin Protectant Prep Pad

- 3: Dust Caps: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp and Fastener
- 1: Catheter Trimmer
- 1: Safety Scalpel: #11
- 1: Chart Sticker
- 1: Checklist/CLIP Sheet
- 1: Flow Rate Information Card with Injection Log
- Patient ID Card
- 1: Patient Information Booklet
- 1: Dressing: STATLOCK®2 Catheter Stabilization Device
- 1: Paper Tape Measure











by ethylene oxide











Not made with natural rubber latex



Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

Contraindications: The Arrow® PICC with Chloragard® Technology is contraindicated for patients with known hypersensitivity to chlorhexidine (refer to product instructions for references).

Arrow International, Inc.
Subsidiary of Teleflex, Incorporated
2400 Bernville Road
Reading, PA 19605 USA



IR-45563-101B (6/16)

¹ Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

² A registered trademark of C. R. Bard, Inc.