# Arterial Catheterization Kit

### **Arrow Select Kits**

## **ASK-00520-AHS**

#### **Contents:**

- 1: Indwelling Catheter: 20 Ga. x 2" (5 cm) Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm)
- (Straight Soft Tip on Both Ends)1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
- 1: Injection Needle: SafetyGlide<sup>™1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- **1:** Syringe: 5 mL Luer-Slip
- **1:** SharpsAway II<sup>™</sup> Locking Disposal Cup
- 1: Drape: 30" x 36" with 3" fenestration, with adhesive
- 1: CSR Wrap
- 1: Needle Holder
- 1: Scissors
- 3: Gauze Pad: 2" x 2"
- 5: Gauze Pad: 4" x 4"
- 1: Dressing: Tegaderm<sup>®2</sup> CHG 8.5 cm x 11.5 cm
- 1: Suture: 3-0 Silk with Curved Needle
- <sup>1</sup> A trademark of Becton, Dickinson and Company.
- <sup>2</sup> A registered trademark of 3M Company.

### > Sharps Safety Features

**INTERNATIONAL** 



Α

#### Rx only.

Select Kits

Arrow

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Contraindications: Inadequate collateral arterial circulation.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc. 2400 Bernville Road, Reading, PA 19605 USA www.arrowintl.com IS-00520-100A (4/15)

Arterial Catheterization Kit







> Sharps Safety Features

ASK-00520-AHS

