



ARROW

# Arterial Catheterization Kit

REF **SAC-01218-X1A**

**18** Ga. | **12** cm catheter length | **.025** inch dia. spring-wire guide

**ErgoPack**

## Contents:

- 1: Indwelling Catheter: 18 Ga. x 12 cm Radiopaque PEBA
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: SafetyGlide™ 22 Ga. x 1-1/2" (3.81 cm) RW
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
- 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)

- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Checklist Sheet
- 1: Medication Label: 1% Lidocaine
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm® 10 cm x 12 cm
- 1: Surgical Apparel: Impervious Gown

- 1: Dressing: BIOPATCH®<sup>4</sup>
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: HemoHopper® Fluid Receptacle

<sup>1</sup>A trademark of Becton, Dickinson and Company.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

<sup>3</sup>A registered trademark of 3M Company.

<sup>4</sup>A registered trademark of Johnson & Johnson Corporation.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

~~LATEX~~

Not made with natural rubber latex



Store between 20-25°C (68-77°F).