



Arterial Catheterization Kit

REF **SAC-01618-1A**

18 Ga. | **16** cm catheter length | **.025** inch dia. spring-wire guide

Contents:

- | | | |
|--|---|---|
| 1: Indwelling Catheter: 18 Ga. x 16 cm Radiopaque PEBA | 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®2} One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution) | 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: SharpsAway [®] Disposal Cup | 1: Dressing: Tegaderm ^{®3} 10 cm x 12 cm |
| 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW | 1: SharpsAway II™ Locking Disposal Cup | 1: Surgical Apparel: Impervious Gown |
| 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm) | 1: Drape: 41" x 55" with 4" fenestration | 1: Surgical Apparel: Bouffant Cap |
| 1: Injection Needle: SafetyGlide™ 22 Ga. x 1-1/2" (3.81 cm) RW | 1: Towel: 13" x 18" (33 cm x 46 cm) | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Needle Holder | 1: Suture: 3-0 Silk with Curved Needle |
| 1: Syringe: 3 mL Luer-Lock | 1: Safety Scalpel: #11 | |
| 1: Syringe: 5 mL Luer-Slip | 1: Medication Label: 1% Lidocaine | ¹ A trademark of Becton, Dickinson and Company. |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution) | 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) | ² A registered trademark of CareFusion or one of its subsidiaries. |
| | | ³ A registered trademark of 3M Company. |

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



Not made with natural rubber latex



Store between 20-25°C (68-77°F).