

REF ASK-04020-MHS

ARROW®



Radial Artery Catheterization Kit

20^{Ga.} | 4.45^{cm} catheter length | .018^{inch dia.} spring-wire guide

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 26 Ga. x 5/8" (1.59 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Skin Protectant Prep Pad (non-sterile solution)
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: Wing Clip: Suture
- 1: SharpsAway II™ Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Safety Scalpel: #11
- 1: Dale® Bendable Armboard[†]
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®[‡] 10 cm x 12 cm
- 1: Suture: 3-0 Braided Silk with Straight Needle

[†]A registered trademark of CareFusion or one of its subsidiaries.

[‡]A registered trademark of Dale Medical Products, Inc.

[§]A registered trademark of 3M Company

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Not made with natural rubber latex



[†]MR Conditional



Store between 20-25°C (68-77°F).

REF ASK-04020-MHS

ARROW®



Radial Artery Catheterization Kit

20^{Ga.} | 4.45^{cm} catheter length | .018^{inch dia.} spring-wire guide

Manufacturer:
Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

IS-04020-115A, Rev. 01 (2020-12)