



Central Venous Catheterization Set



ES-04522-BZ



Conjunto de Cateterização Venosa Central Manterafastado da luz solar. Manterseco.



Conjunto de cateterización venosa centralManteneralejado de la luz del sol. Manteneren un lugarseco.

Contents

- 1: Indwelling Catheter: 22 Ga. x 10 cm Radiopaque Polyurethane with Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 2: Tissue Dilator
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.







Consult instructions for use

Do not reuse

Do not resterilize

STERILE EO

by ethylene oxide







Do not use if package is damaged



Not made with natural rubber latex



Manufacturer Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone Co. Westmeath, Ireland



Catheter Diameter: 22 Ga., 0.9 mm

Priming Volume: 0.20 mL

Priming volumes are approximate and are done without accessories.

Flow Rate: 860 mL/hr Based on 100 cm head height.

IS-04522-131A (3/17)