



# Percutaneous Sheath Introducer Kit for use with 7 - 7.5 Fr. Catheters

REF **ES-09807**

**8** Fr. | **10** cm sheath length | **.035** inch dia. spring-wire guide

**Contents:**

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>1: Sheath: 8 Fr. x 10 cm Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator</li> <li>1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer</li> <li>1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle</li> <li>1: Injection Needle: 25 Ga. x 1" (2.54 cm)</li> <li>1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW</li> <li>1: Local Anesthetic Component</li> <li>1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe</li> <li>1: Syringe: 3 mL Luer-Lock</li> <li>1: Syringe: 5 mL Luer-Slip</li> <li>1: 5 mL 1% Lidocaine HCl Solution</li> </ul> | <ul style="list-style-type: none"> <li>1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)</li> <li>1: Obturator: 8 Fr.</li> <li>1: Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape</li> <li>1: SharpsAway® Disposal Cup</li> <li>1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration, with adhesive</li> <li>2: Towel</li> <li>1: Scalpel: #11</li> <li>5: Gauze Pad: 4" x 4" (10 cm x 10 cm)</li> <li>1: Suture: 1-0 Silk with Straight Needle</li> </ul> |
|---|--|

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

								This Product Contains Dry Natural Rubber.	Store between 20 - 25°C (68 - 77°F).

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

**California Prop. 65**

**Warning: Cancer and Reproductive Harm [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).**

Fluid path components are non-pyrogenic.

DuPont™ Tyvek®  
packaging

**Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IS-09807-100J (03/18)



# Percutaneous Sheath Introducer Kit for use with 7 - 7.5 Fr. Catheters

REF **ES-09807**

**8** Fr. | **10** cm sheath length | **.035** inch dia. spring-wire guide