

# ARROW®



## Central Venous Catheterization Set with Blue FlexTip® Arrowg+ard Blue® Catheter

REF **CS-24706-E**

**14** <sup>Ga.</sup> | **16** <sup>cm</sup> catheter length | **.032** <sup>inch dia.</sup> spring-wire guide

### Contents:

- 1: Indwelling Catheter: 14 Ga. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment<sup>1</sup>, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
- 1: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup>Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.


- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip Arrowg+ard Blue
- de** Besteck für zentralen Venenkatheter mit Blue-FlexTip-Arrowg+ard-Blue-Katheter
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip Arrowg+ard Blue
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip Arrowg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip, Arrowg+ard Blue
- ru** Комплект для катетеризации главных вен с катетером Blue FlexTip Arrowg+ard Blue
- sl** Komplet za centralno venosko katetrizacijo s katetrom Blue FlexTip Arrowg+ard Blue
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip Arrowg+ard Blue
- sv** Set för central venkatetrising med Blue FlexTip Arrowg+ard Blue-kateter
- tr** Santral Venöz Kateterizasyon Seti, Blue FlexTip Arrowg+ard Blue Kateter ile



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.

**Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.48	5400

\* Priming volumes are approximate and are done without accessories.  
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

**Manufacturer:**  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA



IS-24706-130C (06/19)



## Central Venous Catheterization Set with Blue FlexTip® Arrowg+ard Blue® Catheter

REF **CS-24706-E**

**14** <sup>Ga.</sup> | **16** <sup>cm</sup> catheter length | **.032** <sup>inch dia.</sup> spring-wire guide

# ARROW®