

## PICC Seldinger Conversion Set

to be used with 6 Fr. Arrow® PICC

REF **ASK-01663-MST**

**Contents:**

- 1: GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr. Dilator
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW



Caution



Consult instructions for use



Do not reuse



Do not resterilize

STERILE EO

Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.



LOT



DuPont™ Tyvek® packaging

**Teleflex®**

Arrow International, Inc.  
Subsidiary of Teleflex, Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

[www.teleflex.com](http://www.teleflex.com)

IT-01663-100B (4/16)

## PICC Seldinger Conversion Set

REF **ASK-01663-MST**

To be used with **6** Fr. Arrow® PICC

**.018** inch dia. spring-wire guide