



Nitinol Marked Spring-Wire Guide with Arrow® GlideWheel™ Wire Advancer

REF **AW-04432-N**

ARROW

Contents: _____

1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer

Maximum OD per ISO 11070 is 0.82 mm.

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Caution



Consult instructions for use



Do not reuse



Do not sterilize

STERILE EO

Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

LOT



Use By

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com

1W-04432-104A (04/20)



Nitinol Marked Spring-Wire Guide with Arrow® GlideWheel™ Wire Advancer

REF **AW-04432-N**

ARROW

Contents: _____

1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer

Maximum OD per ISO 11070 is 0.82 mm.

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Caution



Consult instructions for use



Do not reuse



Do not sterilize

STERILE EO

Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

LOT



Use By

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com

1W-04432-104A (04/20)