


REF CS-25123-F

# Large-Bore Multi-Lumen Central Venous Catheterization Set

with Blue FlexTip® ARROWg<sup>+</sup>ard Blue® Catheter for High Volume Infusions

# ARROW®



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.40	2457
Medial (12 Ga.)	1.30	10151
Proximal (12 Ga.)	1.20	11039

\* Priming volumes are approximate.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

#### Contents:

- 1:..... Multi-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg<sup>+</sup>ard Antimicrobial Surface Treatment<sup>†</sup>, Extension Line Clamps
- 1:..... Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1:..... Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1:..... Pressure Transduction Probe
- 1:..... Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1:..... Tissue Dilator: Stepped
- 1:..... Drape: 24" x 36"
- 1:..... Scalpel: #11

<sup>†</sup> Licensed under US Patent No. 6,706,024.

#### Rx only.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

**Contraindications:** The ARROWg<sup>+</sup>ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instruction for references).

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Store below 25°C (77°F).

Avoid excessive heat above 40°C (104°F).



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

TEC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
0086 Co. Westmeath, Ireland

**FR** Jeu de cathétérisme veineux central de grand calibre à lumières multiples avec cathéter Blue FlexTip ARROWg<sup>+</sup>ard Blue pour les perfusions de grand volume

Conservez le produit à une température inférieure à 25°C.

Veillez à ne pas exposer le produit à une chaleur excessive supérieure à 40°C.

**DE** Besteck für viellumigen, zentralen, großkalibrigen Venenkatheter mit Blue-FlexTip-ARROWg<sup>+</sup>ard-Blue-Katheter für Infusionen großer Flüssigkeitsmengen

Unter 25°C lagern. Zu hohe Temperaturen über 40°C vermeiden.

**IT** Set per cateterismo della vena centrale multilume di diametro ampio con catetere Blue FlexTip ARROWg<sup>+</sup>ard Blue per volumi di infusione elevati

Conservare ad una temperatura inferiore a 25°C. Evitare il calore eccessivo superiore a 40°C.

**PL** Wielokanalowy zestaw o dużej średnicy do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWg<sup>+</sup>ard Blue do przetaczania dużych objętości płynów infuzyjnych

Przechowywać w temperaturze poniżej 25°C. Chronić przed temperaturą powyżej 40°C.

**PT** Conjunto de Cateterização Venosa Central de Grande Calibre, com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg<sup>+</sup>ard Blue, para Infusões de Grandes Volumes

Conservar a temperaturas inferiores a 25°C. Evitar temperaturas superiores a 40°C.

**RU** Комплект для катетеризации главных вен с многопросветным катетером большого внутреннего диаметра с Blue FlexTip ARROWg<sup>+</sup>ard Blue для вливаний большого объема

Хранить при температуре ниже 25°C. Избегать нагрева свыше 40°C.

**SL** Komplet za centralno venosko katetrizacijo z velikimi luknjami in več svetlinami s katetrom Blue FlexTip ARROWg<sup>+</sup>ard Blue za infundiranje velike prostornine

Shranjujte pod 25°C. Izogibajte se prekomerni vročini nad 40°C.

**ES** Conjunto de cateterización venosa central de luz múltiple y gran diámetro con catéter Blue FlexTip ARROWg<sup>+</sup>ard Blue para infusiones de gran volumen

Almacenar a temperaturas inferiores a 25°C. Evite temperaturas excesivas por encima de 40°C.

**SV** Set för storhålig flerkanalig central venkatetrering med Blue FlexTip ARROWg<sup>+</sup>ard Blue-kateter för infusioner med hög volym

Förvaras vid temperatur under 25°C. Undvik höga temperatur över 40°C.

**TR** Büyük Açıklıklı Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWg<sup>+</sup>ard Blue Kateter ile, Yüksek Hacimli İnfüzyonlar için

25°C'ın altında saklayın. 40°C'ın üzerindeki aşırı ısıdan sakının.

IZ-25123-200A (05/19)

ARROW®  
INTERNATIONAL

REF CS-25123-F 12 Fr. 3 Lumen 20 cm catheter length .035 inch dia. spring-wire guide

Large-Bore Multi-Lumen Central Venous Catheterization Set with ARROWg<sup>+</sup>ard Blue® Catheter

Real to OPEN