



Quad-Lumen Central Venous Catheterization Kit




REF **ASK-42854-PRJ**

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 6" (16 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, ARROWgard Blue PLUS® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: Eclipse™² 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe
- 1: Injection Needle: Eclipse™² 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe
- 1: Syringe: 10 mL Luer-Lock
- 1: Tissue Dilator
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Forceps
- 1: Safety Scalpel: #11
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Flow Rate Information Card with Injection Log
- 2: Prep Sponge Swab
- 2: Gauze Pad: 2" x 2"
- 5: Gauze Pad: 4" x 4"
- 1: Dressing: Tegaderm®³ 10 cm x 12 cm
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tray: Prep
- 1: HemoHopper® Fluid Receptacle

¹ Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.
² A trademark of Becton, Dickinson and Company.
³ A registered trademark of 3M Company.



Lumen	Priming Volume* (mL)	Flow Rate [†] (mL/hr)	MAX Pressure Injection Flow Rate [‡] (mL/sec)
Distal (16 Ga.)	0.38	3350	5
Medial (14 Ga.)	0.45	6000	10
Medial 2 (18 Ga.)	0.35	1750	-
Proximal (18 Ga.)	0.38	1950	-

* Priming volumes are approximate and are done without an injection cap.
[†] Based on 100 cm head height.
[‡] Pressure injection flow rates are determined using media of 11.8 centipoise viscosity, with 152 cm administration set tubing.



Caution Consult instructions for use Do not reuse Do not resterilize Sterilized by ethylene oxide Keep away from sunlight Keep dry Do not use if package is damaged Not made with natural rubber latex Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F)

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. Sterile. Fluid path components are non-pyrogenic.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Contraindications: The ARROWgard Blue Plus® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

J-42854-100A (12/15)



Arrow International, Inc.
 Subsidiary of Teleflex, Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA
www.teleflex.com



DuPont™ Tyvek®
 DuPont™ Tyvek®
 DuPont™ Tyvek®



Quad-Lumen Central Venous Catheterization Kit

REF **ASK-42854-PRJ**

8.5 Fr.

4 Lumen

16 cm catheter length

.032 inch dia. spring-wire guide

