

ARROW



REF CV-12122-F (IPNXXXXXX)

12 Fr. | 2 Lumen | 16 cm catheter length | .035 inch dia. spring-wire guide

## Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

### Contents:

- 1: Two-Lumen Indwelling Catheter: 12 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.


**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

**Contraindications: The Arrow® large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.**

- fr** Jeu de cathétérisme pour hémodialyse à deux lumières avec cathéter Blue FlexTip pour les perfusions de grand volume  
Contre-indications : Le cathéter à deux lumières de gros diamètre Arrow n'est pas conçu pour une hémodialyse à long terme ni pour une utilisation chez des patients souffrant de thromboses vasculaires.
- de** Besteck für zweilumigen Hämodialysekatheter mit Blue FlexTip Katheter für Infusionen großer Flüssigkeitsmengen  
Kontraindikationen: Der doppel- und großlumige Katheter von Arrow ist nicht für die langzeitige Hämodialyse oder zur Verwendung bei Patienten mit thrombosierten Gefäßen geeignet.
- it** Set per cateterismo a due lumi in emodialisi con catetere Blue FlexTip per volumi di infusione elevati  
Controindicazioni - Il catetere Arrow a due lumi di ampio diametro non è progettato per emodialisi a lungo termine o per pazienti che presentano vasi trombizzati.
- pl** Dwukanałowy zestaw cewnika do hemodializy z cewnikiem Blue FlexTip do przetaczania dużych objętości płynów infuzyjnych  
Przeciwwskazania: Dwukanałowy cewnik o dużej średnicy Arrow nie jest przeznaczony do przewlekłej hemodializy ani do stosowania u pacjentów z zakrzepicą naczyń.
- pt** Conjunto de Cateterização para Hemodiálise com Lúmen Duplo e Cateter Blue FlexTip, para Infusões de Grandes Volumes  
Contra-indicações: O cateter de grande calibre e lúmen duplo Arrow não se destina a hemodiálise prolongada nem a utilização em doentes com vasos com trombos.
- ru** Двухпросветный комплект для гемодиализной катетеризации с катетером Blue FlexTip для вливания больших объемов  
Противопоказания. Катетер Arrow с двумя широкими просветами не предназначен для длительного гемодиализа или для использования у пациентов с тромбированными сосудами.
- sl** Komplet za hemodializno katetrizacijo z dvema svetlinama s katetrom Blue FlexTip za infundiranje velike prostornine  
Kontraindikacije: Kateter Arrow z dvema svetlinama in velikimi odprtinami ni zasnovan za dolgoročno hemodializo ali za uporabo pri bolnikih s tromboznimi žilami.
- es** Conjunto de cateterización de dos luces para hemodiálisis con catéter Blue FlexTip para infusiones de gran volumen  
Contraindicaciones: El catéter de doble luz y gran diámetro interior Arrow no ha sido diseñado para la hemodiálisis a largo plazo ni para ser utilizado en pacientes con vasos trombosados.
- sv** Set för tvåkanalig hemodialyskatetrering med Blue FlexTip kateter för infusioner med hög volym  
Kontraindikationer: Arrow grov tvåkanalig kateter är inte avsedd för hemodialys under längre tid eller för användning hos patienter med trombotiserade kärl.
- tr** İki Lümenli Hemodiyaliz Kateterizasyon Seti, Blue FlexTip Kateterli, Yüksek Hacimli İnfüzyonlar için  
Kontrendikasyonlar: Arrow büyük açıklıklı iki lümenli kateter uzun dönemli hemodiyaliz veya trombozlu damarları olan hastalarda kullanım için tasarlanmamıştır.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (12 Ga.)	1.4	220	100
Proximal (12 Ga.)	1.3	370	200
		220	100
		370	200

\* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.



LOT LotNumber

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

CE 0086 Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Teleflex®  
www.teleflex.com

LBL010057 (04/18)



REF CV-12122-F 12 Fr. | 2 Lumen | 16 cm catheter length | .035 inch dia. spring-wire guide ARROW

Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions