

ARROW®



REF CD-13902

(IPNXXXXXX)

9 Fr. | 2 Lumen | 13 cm catheter length | .025 inch dia. spring-wire guide

Pediatric Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

Contents:

- 1: Two-Lumen Indwelling Catheter: 9 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Syringe: 3 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: Scalpel: #11

All components are CE 0086 unless otherwise noted.


Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

Contraindications: The Arrow® large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.

- fr** Jeu de cathétérisme pour hémodialyse pédiatrique à deux lumières avec cathéter Blue FlexTip pour les perfusions de grand volume. Contre-indications : Le cathéter à deux lumières de gros diamètre Arrow n'est pas conçu pour une hémodialyse à long terme ni pour une utilisation chez des patients souffrant de thromboses vasculaires.
- de** Besteck für zweilumigen Hämodialysekatheter für Kinder mit Blue FlexTip Katheter für Infusionen großer Flüssigkeitsmengen. Kontraindikationen: Der doppel- und großlumige Katheter von Arrow ist nicht für die langzeitige Hämodialyse oder zur Verwendung bei Patienten mit thrombosierten Gefäßen geeignet.
- it** Set per cateterismo in emodialisi a due lumi pediatrico con catetere Blue FlexTip per volumi di infusione elevati. Controindicazioni - Il catetere Arrow a due lumi di ampio diametro non è progettato per emodialisi a lungo termine o per pazienti che presentano vasi trombizzati.
- pl** Pediacyjny, dwukanałowy zestaw cewnika do hemodializy z cewnikiem Blue Flex Tip do przetaczania dużych objętości płynów infuzyjnych. Przeciwwskazania: Dwukanałowy cewnik o dużej średnicy Arrow nie jest przeznaczony do przewlekłej hemodializy ani do stosowania u pacjentów z zakrzepicą naczyń.
- pt** Conjunto Pediátrico de Cateterização para Hemodiálise com Lúmen Duplo e Cateter Blue FlexTip, para Infusões de Grandes Volumes. Contra-indicações: O cateter de grande calibre e lúmen duplo Arrow não se destina a hemodiálise prolongada nem a utilização em doentes com vasos com trombos.
- ru** Педиатрический набор для двухпросветной гемодиализной катетеризации с катетером Blue FlexTip для инфузии больших объемов жидкостей. Противопоказания. Катетер Arrow с двумя широкими просветами не предназначен для длительного гемодиализа или для использования у пациентов с тромбированными сосудами.
- sl** Pediatrični komplet za hemodializno katetризacijo z dvema kanaloma in s katetrom Blue FlexTip za infundiranje velike prostornine. Kontraindikacije: Kateter Arrow z dvema svetlinama in velikimi odprtinama ni zasnovan za dolgoročno hemodializo ali za uporabo pri bolnikih s tromboznimi žilami.
- es** Conjunto de cateterización de dos luces para hemodiálisis para uso pediátrico con catéter Blue FlexTip para infusiones de gran volumen. Contraindicaciones: El catéter de doble luz y gran diámetro interior Arrow no ha sido diseñado para la hemodiálisis a largo plazo ni para ser utilizado en pacientes con vasos trombosados.
- sv** Set för pediatrik tvåkanalig hemodialyskatetrisering med Blue FlexTip kateter för infusioner med hög volym. Kontraindikationer: Arrow grov tvåkanalig kateter är inte avsedd för hemodialys under längre tid eller för användning hos patienter med trombotiserade kärl.
- tr** Pediatrik İki Lümenli Hemodializ Kateterizasyon Seti, Blue FlexTip Kateterli, Yüksek Hacimli İnfüzyonlar İçin. Kontrendikasyonlar: Arrow büyük açıklıklı iki lümenli kateter uzun dönemli hemodializ veya trombozlu damarları olan hastalarda kullanım için tasarlanmamıştır.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (14 Ga.)	0.64	84 189	100 200
Proximal (14 Ga.)	0.55	84 189	100 200

* Priming volumes are approximate and are done without accessories.
 † Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.

Lot Number
 YYYY-MM-DD

(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (10)LotNumber

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA
Teleflex®
 www.teleflex.com
 LBL010063 (04/18)



REF CD-13902

9 Fr. | 2 Lumen | 13 cm catheter length | .025 inch dia. spring-wire guide

Pediatric Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

ARROW®