



REF CV-15122-UF

(IPNXXXXXX)

12 Fr. | 2 Lumen | 20 cm catheter length | .035 inch dia. spring-wire guide

You-Bend™ Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

Contents:

- 1: Two-Lumen Indwelling Catheter‡: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.


‡MR Conditional.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

Contraindications: The Arrow® large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.

- fr** Jeu de cathétérisme pour hémodialyse à deux lumières You-Bend avec cathéter Blue FlexTip pour les perfusions de grand volume. Contre-indications : Le cathéter à deux lumières de gros diamètre Arrow n'est pas conçu pour une hémodialyse à long terme ni pour une utilisation chez des patients souffrant de thromboses vasculaires.
- de** Besteck für zweilumigen You-Bend-Hämodialysekatheter mit Blue FlexTip Katheter für Infusionen großer Flüssigkeitsmengen. Kontraindikationen: Der doppel- und großlumige Katheter von Arrow ist nicht für die langzeitige Hämodialyse oder zur Verwendung bei Patienten mit thrombosierten Gefäßen geeignet.
- it** Set per cateterismo in emodialisi a due lumi You-Bend con catetere Blue FlexTip per volumi di infusione elevati. Controindicazioni - Il catetere Arrow a due lumi di ampio diametro non è progettato per emodialisi a lungo termine o per pazienti che presentano vasi trombizzati.
- pt** Conjunto de Cateterização para Hemodiálise You-Bend com Lúmen Duplo e Cateter Blue FlexTip, para Infusões de Grandes Volumes. Contra-indicações: O cateter de grande calibre e lúmen duplo Arrow não se destina a hemodiálise prolongada nem a utilização em doentes com vasos com trombos.
- es** Conjunto de cateterización de dos luces para hemodiálisis You-Bend con catéter Blue FlexTip para infusiones de gran volumen. Contraindicaciones: El catéter de doble luz y gran diámetro interior Arrow no ha sido diseñado para la hemodiálisis a largo plazo ni para ser utilizado en pacientes con vasos trombosados.
- sv** You-Bend set för tvåkanalig hemodialyskaterisering med Blue FlexTip kateter för infusioner med hög volym. Kontraindikationer: Arrow grov tvåkanalig kateter är inte avsedd för hemodialys under längre tid eller för användning hos patienter med trombotiserade kärl.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (12 Ga.)	1.8	200	100
Proximal (12 Ga.)	1.7	340	200
		200	100
		340	200

* Priming volumes are approximate and are done without accessories.
 † Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.



LOT LotNumber
 YYY-YY-MM-DD



(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (10)LotNumber

ECI Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Teleflex®
 www.teleflex.com
 LBL010064 (04/18)