


**REF SA-15122-F**

(IPNXXXXXX)

**12** Fr.**2** Lumen**20** cm catheter length**.035** inch dia. spring-wire guide**ARROW®****Two-Lumen Hemodialysis Catheterization Kit with Blue FlexTip® Catheter for High Volume Infusions****Contents:**

- 1: Two-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 1: Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive
- 2: Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive
- 1: Safety Scalpel: #11
- 2: Prep Sponge Swab
- 2: Valve: MicroCLAVE®\* Neutral Displacement Connector

\*A registered trademark of ICU Medical, Inc.

**All components are CE 0086 unless otherwise noted.**

Catheter Cross Section

Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (12 Ga.)	1.6	200 340	100 200
Proximal (12 Ga.)	1.5	200 340	100 200

\* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.** Not made with natural rubber latex. Fluid path components are non-pyrogenic.**Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.****EC REP****CE**  
0086Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, IrelandArrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA**LOT**

LotNumber



YYYY-MM-DD



(01)XXXXXXXXXXXXXX

(17)YYMMDD

(10)LotNumber

**Teleflex®**  
www.teleflex.com

LBL010083 (08/19)

**REF SA-15122-F****LOT** LotNumber

YYYY-MM-DD

(01)XXXXXXXXXXXXXX

**REF SA-15122-F****LOT** LotNumber

YYYY-MM-DD

(01)XXXXXXXXXXXXXX