

REF SA-12802 (IPNXXXXXXX)



6 catheter length

spring-wire guide

## Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

- Two-Lumen Indwelling Catheter: 8 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other)
- with Arrow Advancer
- with Arrow Advancer

  Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
  Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW

  Pressure Transduction Probe
  Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL

  Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm

  SecondSite™ Adjustable Hub: Fastener

  SecondSite™ Adjustable Hub: Catheter Clamp

  Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive

  Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive


  Safety Scalpel: #11

- Safety Scalpel: #11
- Prep Sponge Swab Valve: MicroCLAVE®¹ Neutral Displacement Connector
- <sup>1</sup>A registered trademark of ICU Medical, Inc.

All components are CE 0086 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.



Catheter Cross Section

ļ
[
[ F



(mL/hr)

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.
























(01)XXXXXXXXXXXXXXX (17)YYMMDD (10)LotNumber



Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA







**REF SA-12802** 





