

REF SA-12854 (IPNXXXXXXX)



6 catheter length

spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® **Catheter**

- Quad-Lumen Indwelling Catheter: $8.5 \, \text{Fr.} \times 16 \, \text{cm}$ Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps Spring-Wire Guide, Nitinol, Marked: .032" ($0.81 \, \text{mm}$) dia. $\times 17-3/4"$ ($45 \, \text{cm}$) (Straight Soft Tip on One End "J" Tip on Other)
- with Arrow Advancer
- with Arrow Advancer

 Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
 Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW

 Pressure Transduction Probe
 Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL

 Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm

 SecondSite™ Adjustable Hub: Fastener

 SecondSite™ Adjustable Hub: Catheter Clamp

 Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive

 Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive


 Safety Scalpel: #11

- Safety Scalpel: #11
- Prep Sponge Swab Valve: MicroCLAVE®¹ Neutral Displacement Connector
- ¹A registered trademark of ICU Medical, Inc.

All components are CE 0086 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.



_		
	Catheter	

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.39	3200	
Medial 1 (14 Ga.)	0.47	5400	
Medial 2 (18 Ga.)	0.36	1600	
Proximal (18 Ga.)	0.38	1800	

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.






















YYYY-MM-DD



(01)XXXXXXXXXXXXXXX (17)YYMMDD (10)LotNumber



Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA







REF SA-12854

8.5 Fr.



