



20 cm catheter length

1032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

- Multi-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End "J" Tip on Other)
- with Arrow Advancer with ECG Mark


- with Arrow Advancer with ECG Mark
 Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
 Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
 Pressure Transduction Probe
 Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
 Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
 SecondSite™ Adjustable Hub: Fastener
 SecondSite™ Adjustable Hub: Catheter Clamp
 Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive
 Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive
 Safety Scaloel: #11
- Safety Scalpel: #11
- Prep Sponge Swab Valve: MicroCLAVE®¹ Neutral Displacement Connector
- ¹A registered trademark of ICU Medical, Inc.

All components are CE 0086 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.47	3530	
Medial (14 Ga.)	0.89	7555	
Proximal (16 Ga.)	0.47	3290	

^{*} Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.






















(01)XXXXXXXXXXXXXXX (17)YYMMDD (10)LotNumber



Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA







SA-15853

8.5 Fr.



