





Radial Artery Catheterization Kit

cm catheter length

Contents:

- Catheter: 18 Ga. x 4.45 cm Radiopaque FEP over 20 Ga. TW Introducer Needle with Integral .025" (0.64 mm) dia. Spring-Wire Guide Injection Needle: 25 Ga. x 1" (2.54 cm)

- Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
 Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
 Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
 Syringe: 3 mL Luer-Lock
 Foil Pack 70% Alcohol Prep Pad (non-sterile solution)
 Swabstick Foil Pack Skin Protectant (non-sterile solution)
 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution (non-sterile solution)
- sterile solution)
- Wing Clip: Suture Towel: 13" x 18" (33 cm x 46 cm) Gauze Pad: 2" x 2" (5 cm x 5 cm)
- ¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex

Store between 20 – 25°C (68 – 77°F).

Fluid path components are non-pyrogenic.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.







Radial Artery Catheterization Kit

Contents:

- Catheter: 18 Ga. x 4.45 cm Radiopaque FEP over 20 Ga. TW Introducer Needle with Integral .025" (0.64 mm) dia. Spring-Wire Guide Injection Needle: 25 Ga. x 1" (2.54 cm)

- Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
 Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
 Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
 Syringe: 3 mL Luer-Lock
 Foil Pack 70% Alcohol Prep Pad (non-sterile solution)
 Swabstick Foil Pack Skin Protectant (non-sterile solution)
 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution (non-sterile solution) sterile solution)
- Wing Clip: Suture Towel: 13" x 18" (33 cm x 46 cm) Gauze Pad: 2" x 2" (5 cm x 5 cm)
- ¹A registered trademark of CareFusion or one of its subsidiaries

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex

Store between 20 – 25°C (68 – 77°F). Fluid path components are non-pyrogenic.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.





Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA























Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA





REF AK-04018

Radial Artery Catheterization Kit

4.45 cm catheter length





4.45 cm catheter length

