

ARROW

ARROW

REF AK-04020-SPC

REF AK-04020-SPC



Radial Artery Catheterization Kit with Integral Needle Protection

20 Ga. **4.45** cm catheter length



Radial Artery Catheterization Kit with Integral Needle Protection

20 Ga. **4.45** cm catheter length

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Swabstick Foil Pack Skin Protectant (non-sterile solution)
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Frepp® Solution (non-sterile solution)
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

*A registered trademark of CareFusion or one of its subsidiaries.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Fluid path components are non-pyrogenic.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Swabstick Foil Pack Skin Protectant (non-sterile solution)
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Frepp® Solution (non-sterile solution)
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

*A registered trademark of CareFusion or one of its subsidiaries.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Fluid path components are non-pyrogenic.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.



20°C (68°F) to 25°C (77°F)



20°C (68°F) to 25°C (77°F)



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
LBL015384 (01/19)



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
LBL015384 (01/19)



REF AK-04020-SPC

Radial Artery Catheterization Kit with Integral Needle Protection

20 Ga. **4.45** cm catheter length

ARROW



REF AK-04020-SPC

Radial Artery Catheterization Kit with Integral Needle Protection

20 Ga. **4.45** cm catheter length

ARROW