

ARROW

REF AK-04120



Radial Artery Catheterization

20 Ga. 3.81 cm catheter length

Contents:

- 1: Catheter: 20 Ga. x 3.81 cm Radiopaque Polyurethane over Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad
- 1: Skin Protectant Prep
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Frepp® Solution
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

ARROW

REF AK-04120



Radial Artery Catheterization

20 Ga. 3.81 cm catheter length

Contents:

- 1: Catheter: 20 Ga. x 3.81 cm Radiopaque Polyurethane over Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad
- 1: Skin Protectant Prep
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Frepp® Solution
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
LBL015386 R01 (2022-08)

LOT



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
LBL015386 R01 (2022-08)



REF AK-04120

Radial Artery Catheterization

20 Ga. 3.81 cm catheter length

ARROW



REF AK-04120

Radial Artery Catheterization

20 Ga. 3.81 cm catheter length

ARROW