

# ARROW

REF AK-04122



## Radial Artery Catheterization Kit

22 Ga. 3.49 cm catheter length

### Contents:

- 1: Spring-Wire Guide/Tube Assembly: .015" (0.38 mm) dia. x 4-7/8" (12.4 cm)
- 1: Catheter: 22 Ga. x 1-3/8" (3.49 cm) Radiopaque over 23 Ga. TW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad (non-sterile solution)
- 1: Swabstick Foil Pack Skin Protectant (non-sterile solution)
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Frepp® Solution (non-sterile solution)
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

\*A registered trademark of CareFusion or one of its subsidiaries.

### Rx only.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

### California Prop. 65

**WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Fluid path components are non-pyrogenic.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

# ARROW

REF AK-04122



## Radial Artery Catheterization Kit

22 Ga. 3.49 cm catheter length

### Contents:

- 1: Spring-Wire Guide/Tube Assembly: .015" (0.38 mm) dia. x 4-7/8" (12.4 cm)
- 1: Catheter: 22 Ga. x 1-3/8" (3.49 cm) Radiopaque over 23 Ga. TW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad (non-sterile solution)
- 1: Swabstick Foil Pack Skin Protectant (non-sterile solution)
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Frepp® Solution (non-sterile solution)
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

\*A registered trademark of CareFusion or one of its subsidiaries.

### Rx only.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

### California Prop. 65

**WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Fluid path components are non-pyrogenic.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)  
LBL015388 (01/19)



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)  
LBL015388 (01/19)



## Radial Artery Catheterization Kit

22 Ga. 3.49 cm catheter length

# ARROW



## Radial Artery Catheterization Kit

22 Ga. 3.49 cm catheter length

# ARROW