

# ARROW

REF AK-04020-C



## Radial Artery Catheterization Kit

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length

### Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 26 Ga. x 5/8" (1.59 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad (non-sterile solution)
- 1: Swabstick Foil Pack Skin Protectant (non-sterile solution)
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Frepp® Solution (non-sterile solution)
- 1: Wing Clip: Suture
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

\*A registered trademark of CareFusion or one of its subsidiaries.

### Rx only.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

### California Prop. 65

**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Fluid path components are non-pyrogenic.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

# ARROW

REF AK-04020-C



## Radial Artery Catheterization Kit

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length

### Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 26 Ga. x 5/8" (1.59 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad (non-sterile solution)
- 1: Swabstick Foil Pack Skin Protectant (non-sterile solution)
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Frepp® Solution (non-sterile solution)
- 1: Wing Clip: Suture
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

\*A registered trademark of CareFusion or one of its subsidiaries.

### Rx only.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

### California Prop. 65

**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Fluid path components are non-pyrogenic.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.



20°C (68°F) to 25°C (77°F)

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)  
LBL015390 (01/19)



20°C (68°F) to 25°C (77°F)

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)  
LBL015390 (01/19)



REF AK-04020-C



## Radial Artery Catheterization Kit

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length

# ARROW

REF AK-04020-C



## Radial Artery Catheterization Kit

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length

# ARROW