

ARROW

REF AK-04220



QuickFlash® Radial Artery Catheterization

20 Ga. 3.81 cm catheter length

Contents:

- 1: Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 27 Ga. x 1/2" (1.27 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad
- 1: Skin Protectant Prep
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

ARROW

REF AK-04220



QuickFlash® Radial Artery Catheterization

20 Ga. 3.81 cm catheter length

Contents:

- 1: Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 27 Ga. x 1/2" (1.27 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad
- 1: Skin Protectant Prep
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



20°C (68°F) to 25°C (77°F)



20°C (68°F) to 25°C (77°F)

LOT



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
LBL015392 R02 (2023-03)

LOT



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
LBL015392 R02 (2023-03)

REF AK-04220



QuickFlash® Radial Artery Catheterization

20 Ga. 3.81 cm catheter length

ARROW

REF AK-04220



QuickFlash® Radial Artery Catheterization

20 Ga. 3.81 cm catheter length

ARROW