



## REF MTO-09801-CHA

.035 spring-wire guide **10**<sup>cm</sup> sheath length 8.5<sup>Fr.</sup>

## Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

## Contents:

- Sheath: 8.5 Fr. x 4-1/8" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral 1: Hemostasis Valve/Side Port and Tissue Dilator
- Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on 1:
- One End "J" Tip on Other) with Arrow Advancer Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1:
- Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) Injection Needle: 25 Ga. x 1" (2.54 cm) 1: 1.
- Pressure Transduction Probe 1:
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: Obturator: 8 Fr
- Cath-Gard<sup>®</sup> Catheter Contamination Shield: 80 cm with Sterile Tape 1:
- 1: Dust Cap: Non-Vented
- 1. SharpsAway® Disposal Cup
- Drape: 24" x 36" (61 cm x 91 cm) 1:
- 1: Drape: 24" x 36" with 4" fenestration
- 1. Scalpel: #11
- 1: Stopcock: 3-way High-Flow

## Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

A TI & STERILE EO 🔆 🔶

This product contains dry natural rubber.

Fluid path components are non-pyrogenic.

- 2. Prep Sponge Swab 2:
  - Gauze Pad: 2" x 2" (5 cm x 5 cm) Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 5: 1: Suture: 3-0 Nylon with Straight Needle
- 1: Tray: Prep

Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA





LOT

Ω







Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

LATEX